



## Tattoo and mark-free radiation therapy with enhanced safety, accuracy, and patient comfort

We offer radiation therapy for all cancer types and specialize in Breast Deep Inspiration Breath Hold (DIBH), Brain Stereotactic Radiosurgery (SRS) and Stereotactic Body Radiation Therapy (SBRT).

Our advanced technology, AlignRT®, helps deliver radiation with sub-millimeter accuracy to help protect surrounding healthy tissue and organs.

This latest technology also allows us to offer tattoo and mark-free treatment so that patients can choose to not have permanent reminders.

### Our cutting-edge technology offers significant clinical benefits

AlignRT is a contact-free technology that uses 3D cameras and computer vision to help ensure treatment is delivered accurately, safely and without tattoos and marks. Implementing radiation treatment with AlignRT helps ensure treatment is delivered exactly as planned.

#### Exceptional accuracy

AlignRT tracks patients in real-time during the entire treatment delivery procedure with better than 1mm accuracy.<sup>1,2,3</sup> Compared to traditional methods, AlignRT® provides more setup accuracy.<sup>5,6,7</sup>

#### Automated safety monitoring system

AlignRT automatically pauses the beam if clinically significant movement is detected, or if a patient isn't at the right breath-hold level.<sup>4,5,6,7</sup>

#### Enhanced patient comfort


With AlignRT, there is no need to use immobilization devices or tattoos and marks. AlignRT is clinically proven to be as accurate as solutions relying on additional internal imaging.<sup>8,9</sup>

For Stereotactic Radiotherapy (SRS), AlignRT is a frameless, open-faced mask solution that offers patients enhanced safety and comfort. This may be preferred by patients with claustrophobia.



### PROTECT THE HEART Cardiac perfusion defects at 6 months


As measured using SPECT imaging



### Patients will spend more time and money for a tattooless and markless option

78% of patients would choose treatment which avoided tattoos and/or marks, even if additional efforts required, additional costs, distance, and time for travel.<sup>8</sup>

#### Additional travel, time & cost



## (INSERT CENTER INFORMATION)

### Offering enhanced safety, accuracy, and tattoo and mark-free treatment

We appreciate your referrals.  
Your patients will appreciate our treatment.

- Wen et al. Technical Note: Evaluation of the systematic accuracy of a frameless, multiple image modality guided, linear accelerator based stereotactic radiosurgery system. *Med Phys.* 2016 May;43(5):2527
- Covington E, et al. Submillimeter monitoring of intrafraction patient movement with optical surface imaging. *AAPM Annual Meeting* 2018.
- Vision RT. Advanced Camera Optimization [White paper] 2018
- Leong B et al. Impact of Use of Optical Surface Imaging on Initial Patient Setup for Stereotactic Body Radiotherapy Treatments. *Journal of Applied Clinical Medical Physics* 2019; 20 (12):149-158
- Herron E, et al. Surface guided radiation therapy as a replacement for patient marks in treatment of breast cancer. *International Journal of Radiation Oncology • Biology • Physics* 2018;102 (3):e492-e493
- Allen and Ramtohl. Markerless Radiotherapy at the Queen Elizabeth Hospital Birmingham. *Radiography.* 2020; 26(1):S8-9
- Sarkar V, et al. An evaluation of the consistency of shifts reported by three different systems for non-coplanar treatments. *Journal of radiosurgery and SBRT* 2018;5 (4):323-330.
- Wiant D, et al. Direct comparison between surface imaging and orthogonal radiographic imaging for SRS localization in phantom. *Journal of Applied Clinical Medical Physics.* 2019; 20 (1): 137-144
- Menyhart G. SGRT for Breath Hold SBRT: A Measured Response. *SGRT Annual Meeting, Atlanta* 2016
- Marks et al. The incidence and functional consequences of RT-associated cardiac perfusion defects. *Int J Radiat Oncol Biol Phys.* 2005 Sep 1;63(1):214-23
- Zagar et al. Utility of Deep inspiration breath-hold for left sided breast radiation therapy in preventing early cardiac perfusion defects – A Prospective Study. *Int J Radiat Oncol Biol Phys.* 2017